



THE PERFECTION OF  
NATURAL SOUND

WIDEX **MOMENT SHEER™**

**WIDEX**

SOUND LIKE NO OTHER



# THE PERFECTION OF NATURAL SOUND IS WHEN I FEEL CONNECTED TO MY LIFE AGAIN

We know that the first experience of a hearing aid is vital in setting the direction of the journey to better hearing. A negative experience can leave you wondering if you were right to try a hearing aid at all, but a positive one can make you feel more connected to your life than you have been for years.

We believe that it's possible to give you a sound so natural that you will be simply amazed from the very first moment. How? With Widex PureSound™.

Typically, sound processed in a hearing aid reaches the eardrum later than sound heard directly. And when these two 'out of sync' signals mix, you get an artificial sound. But with PureSound, you can hear every moment like you used to.

## NATURAL SOUND AT ITS BEST

We asked real hearing aid users to wear our hearing aids in their everyday lives and tell us what they thought. The results were overwhelmingly positive, with 95% perceiving the sound as 'natural' and 'clear'.







## STREAM FROM iOS AND ANDROID

All our direct streaming hearing aids can be connected via 2.4 GHz Bluetooth to both iOS and ASHA compatible Android devices, Widex TV Play, Widex Sound Assist, and other DEX accessories.

So, now you're prepared for all the streaming entertainment the world has to offer in high-quality sound.

### WIDEX MOMENT SHEER



sRIC R D

### WIDEX MOMENT Direct Streaming



RIC 312 D



BTE R D



BTE 13 D





## AVAILABLE CHARGERS

### STANDARD CHARGERS

(sRIC R D, BTE R D)

Regardless of which rechargeable Moment hearing aids you use, there is a standard charger available.

- 4-hour charge for full-day use
- 30-min fast charge for 4-hour use
- Slim and light
- LED light indicators to show charging status



### CHARGE N CLEAN CHARGER

- 4-hour charge for full-day use
- 30-min charge for 4-hour HA use
- Sanitizing function of ear ware
- Drying function






## CONTROL AT YOUR FINGERTIPS




With the Widex Moment hearing aids, you have the option to fully control them using the Widex Moment app on your Apple or Android smartphone\*. It has never been easier to fine-tune your hearing, so you get exactly the sound you prefer – in exactly the moments that are important to you.

Rigorous testing has verified that the app is easy and intuitive. And a 'dark mode' ensures discreet use.

\* The Widex Moment app works with the following Widex Moment models: sRIC R D (MRR4D), mRIC R D (MRR2D), RIC 312 D (MRB2D), BTE R D (MBR3D) and BTE 13 D (MBB3D). Use the Widex Tonelink app for other models for control of key functions.



WIDEX **MOMENT** App



WIDEX **TONELINK** App



## THE PERFECTION OF NATURAL SOUND

What if a hearing aid could give you access to the listening experiences of thousands of real users just like you, so they could give you a personally recommended sound in an instant?

With MySound, which is part of the MOMENT App, you can. MySound combines real human experiences with the immense power of artificial intelligence (AI).

### **More control. More insights**

With MySound 2.0, we've added compression to the list of parameters that you can control. This means that you can now create programs that adjust both equalizer and compression in a way that suits your need and moment – delivering an even more truly personalized hearing experience.

## IS WHEN I FEEL LIKE I'M IN CONTROL

### **Made for You**

Simply by selecting the most relevant activity and listening intent, our AI instantly presents you with a recommended sound profile based on the way thousands of other users preferred to hear in a similar situation.

### **Made by You — now even faster**

The AI system guides you to a better sound through a series of comparisons. This allows you to tailor the sound to your exact taste.

Because our calculations now run in the cloud, the 'Made By You' process is now faster than ever. This means you have access to high-speed computations, whether you have the newest smartphone or an older model.

## WIDEX MY GUIDE – YOUR GUIDE TO MAKING THE MOST OF MOMENT

Introducing: My Guide. A section of the Moment App that gives you advice in the form of step-by-step instructions and videos to help you with everything from turning your hearing aids on to personalizing your sound with MySound.

There's also a series of short courses designed to help you get used to daily life with hearing aids.



# WIDEX **MOMENT** MODELS AND STYLES

Widex Moment is available as RIC, BTE and In-the-ear models. Your hearing care professional can help you choose the ideal model for your specific hearing loss and lifestyle.

## Choose your colors for RIC and BTE



## RIC



## BTE



## IN-THE-EAR models



## GET MORE FROM YOUR WIDEX **MOMENT**


























### Wireless accessories – connect, control and communicate

DEX accessories can connect effortlessly to all Moment devices and offer a range of additional functionalities.

TV Play, SoundConnect and Widex Sound Assist are accessories for our Direct streaming devices. TV Play offers TV sound streaming, Sound Connect allows audio streaming from PC/laptops and Widex Sound Assist offers a whole range of functionalities.

# WHY WIDEX **MOMENT** IS RIGHT FOR YOU

To help you hear your best, Widex Moment features exclusive Sound Class Technology that automatically adapts your hearing aid to your specific listening requirements. Without you pushing any buttons, Widex Moment ensures optimal clarity, comfort and audibility in each different situation.

PERFORMANCE LEVEL/BENEFIT	SOCIAL	MUSIC		QUIET		PARTY		TRANSPORT		URBAN	
		-	Classical or contemporary	-	with speech	-	with speech	-	with speech	-	with speech
Widex Moment 440											
Widex Moment 330											
Widex Moment 220											
Widex Moment 110											

PERFORMANCE LEVEL/BENEFIT	Improved speech understanding Get accentuated speech and reduced noise from the surroundings	Wind noise reduction Hear clearly even when it's windy	Natural sound experience Choose between the renowned Widex sound or the new Widex PureSound	Automatic sound adjustment No need to adjust the hearing aid in every situation	Direct Bluetooth® streaming* Stream sound directly from Apple or Android devices	Easy to control* Intuitive app and artificial Intelligence guide you to the best sound	Widex Remote Care Enables fitting and fine-tuning at any location the client requests	Real-life insights Facilitates counselling by providing access to the clients' personal programs
Widex Moment 440	● ● ● ● ● ●	● ● ● ● ● ●	● ● ● ● ● ●	● ● ● ● ● ●	YES	YES	YES	YES
Widex Moment 330	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	YES	YES	YES	YES
Widex Moment 220	● ● ● ●	● ● ● ●	● ● ● ●	● ● ● ●	YES	YES	YES	YES
Widex Moment 110	● ●	● ●	● ●	● ●	YES **	YES **	YES	YES

\* mRIC R D, sRIC R D, RIC 312 D, BTE 13 D, BTE R D  
\*\* Doesn't include MySound





SOUND LIKE NO OTHER



Widex is a registered trademark of WIDEX A/S.  
Printed by Rosendahls a/s. 2022-11. 9 502 5335 001 #06

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.  
Android, Google Play and Google Play logo are trademarks of Google Inc.  
Subject to change with no further notice.

[WWW.WIDEX.COM](http://WWW.WIDEX.COM)